

Product Datasheet

tuba8l2 antibody (Biotin) (orb53158)

Description	tuba8l2 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish Tubulin alpha chain protein (1-449AA)
Target	tuba8l2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6PC95
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.