

Product Datasheet

glcB antibody (HRP) (orb53220)

Description	glcB antibody (HRP)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli Malate synthase G protein (2-723AA)
Target	glcB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P37330
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.