

Product Datasheet

Propargyl-choline (orb532316)

Description	Propargyl-choline
Form/Appearance	solid or clear oil; Color: white to off-white
Storage	store at -20 °C
Note	For research use only
Application notes	Choline-containing phospholipid synthesis monitoring in cell culture and whole organisms.
Formula	C ₇ H ₁₄ NO
Purity	H NMR conforms to structure
MW	Theoretical MW: 128.19 g/mol (cation)
SMILES	<chem>OCC[N+](C)(CC#C)C[NO NAME]</chem>
Solubility (25°C)	water, 1x PBS buffer
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.