

Product Datasheet

WSS1 antibody (FITC) (orb53237)

Description	WSS1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	FITC
Immunogen	Recombinant Saccharomyces cerevisiae DNA-dependent metalloprotease WSS1 protein (1-269AA)
Target	WSS1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P38838
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.