


Product Datasheet

OSM34 antibody (orb53247)

Description	OSM34 antibody
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Arabidopsis thaliana Osmotin-like protein OSM34 protein (23-244AA)
Target	OSM34
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P50700
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Recombinant pro...