

Product Datasheet

T4 DNA Ligase (orb532701)

Description	T4 DNA Ligase catalyzes the formation of a phosphodiester bond between juxtaposed 5' phosphate and 3...
Form/Appearance	liquid (Supplied in 10 mM Tris-HCl pH 7.4, 50 mM KCl, 0.1 mM EDTA, 1 mM DTT, 200 µg/ml rAlbumin and 50% [v/v] glycerol)
Concentration	2.5 Weiss units/µl (500 CE units/µl)
Storage	store at -20 °C
Note	For research use only
Expiration Date	6 months from date of receipt.