

Product Datasheet

Anionic trypsin antibody (Biotin) (orb53286)

Description	Anionic trypsin antibody (Biotin)
Species/Host	Rabbit
Reactivity	Canine
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Dog Anionic trypsin protein (24-247AA)
Target	Anionic trypsin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P06872
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.