

Product Datasheet

Low Range DNA Ladder (orb532916)

Description	Low Range DNA Ladder consists of 11 DNA fragments ranging from 50 bp to 1 kb. It is designed to show...
Form/Appearance	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 8% (w/v) Glycerol, Tartrazine and Xylene Cyanol FF; Color: green
Concentration	120 ng/μl
Storage	store at -20 °C
Note	For research use only
Application notes	Recommended Load: 5 μl per lane.
Expiration Date	12 months from date of receipt.