

Product Datasheet

50 bp DNA Ladder (orb532921)

Description	50 bp DNA Ladder consists of 10 DNA fragments ranging from 50 bp to 500 bp and is supplied in ready-...
Form/Appearance	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 10% (w/v) Glycerol and Orange G; Color: orange
Concentration	105 ng/μl
Storage	store at -20 °C
Note	For research use only
Application notes	Recommended Load: 5 μl per lane.
Expiration Date	12 months from date of receipt.