

Product Datasheet

500 bp DNA Ladder (orb532924)

Description	500 bp DNA Ladder consists of 10 DNA fragments ranging from 500 bp to 5 kb and is supplied in ready-...
Form/Appearance	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 10% (w/v) Glycerol, Bromophenol Blue and Xylene Cyanol FF; Color: blue
Concentration	105 ng/μl
Storage	store at -20 °C
Note	For research use only
Application notes	Recommended Load: 5 μl per lane.
Expiration Date	12 months from date of receipt.