

Product Datasheet

Fluorescent High Range DNA Ladder (orb532928)

Description	Fluorescent High Range DNA Ladder consists of 12 DNA fragments ranging from 500 bp to 10 kb. It is d...
Form/Appearance	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 10% (w/v) Glycerol, SYBR® Green, Bromophenol Blue and Xylene Cyanol FF; Color: blue
Concentration	130 ng/μl
Storage	store at -20 °C
Note	For research use only
Application notes	Recommended Load: 5 μl per lane.
Expiration Date	12 months from date of receipt.