

Product Datasheet

Fluorescent 100 bp DNA Ladder (orb532932)

Description	Fluorescent 100 bp DNA Ladder consist of 10 DNA fragments ranging from 100 bp to 1 kb. The prestaine...
Form/Appearance	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 10% (w/v) Glycerol, SYBR® Tartrazine and Xylene Cyanol FF; Color: green
Concentration	105 ng/μl
Storage	store at -20 °C
Note	For research use only
Application notes	Recommended Load: 5 μl per lane.
Expiration Date	12 months from date of receipt.