

Product Datasheet

Low MWCO Filtration Plate (orb533000)

Description Plate Material: Polypropylene, hydrophobic Number of Wells: 96 Well Capacity: 400 µl Filter Membran...


Storage store at ambient temperature

Note For research use only

Application notes The plates are used to separate proteins/protein detergent complexes by their sizes (diameters). Polyethersulfone membranes have very low non-specific binding to proteins and thus optimize yields and purity. Well Capacity: 400 µl Filter Membrane: Polyethersulfone (PES) Molecular Weight Cutoff: 100 KDa.

Hazard Information Non-Toxic

Expiration Date 12 months from date of receipt.



Low MWCO
Filtration
Plate.