

Product Datasheet

N1-Methyl-AMP (orb533054)

Description	N1-Methyl-AMP
Form/Appearance	solid; Color: white to off-white
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{max} 257 nm, ϵ 14.6 L mmol ⁻¹ cm ⁻¹ (Tris-HCl, pH 7.5).
Formula	C ₁₁ H ₁₆ N ₅ O ₇ P
Purity	≥ 95% (HPLC)
MW	Theoretical MW: 361.25 g/mol (free acid); Detected MW: 361.08 g/mol (free acid)
SMILES	<chem>O[C@H]1[C@H](n2cnc3c(=N)n(cnc23)C)O[C@H](COP(=O)(O)O)[C@H]1O</chem>
Expiration Date	12 months from date of receipt.