

Product Datasheet

IppNHp (orb533085)

Description	IppNHp
Form/Appearance	solid; Color: colorless to white
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{max} 249 nm, ϵ 12.2 L mmol ⁻¹ cm ⁻¹ (Tris-HCl, pH 7.0).
Formula	C ₁₀ H ₁₆ N ₅ O ₁₃ P ₃
Purity	≥ 95% (HPLC)
MW	Theoretical MW: 507.18 g/mol (free acid); Detected MW: 507.00 g/mol (free acid)
SMILES	[C@@H]1([C@H]([C@@H](O[C@@H]1COP(=O)(OP(=O)(NP(=O)(O)O)O)O)n1c2c(nc1)c(=O)[nH]cn2)O)O
Expiration Date	12 months from date of receipt.