

## Product Datasheet

### Crystal Hot Start Master (2x) (orb533330)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Crystal Hot Start Master is a 2 x conc. ready-to-use master mix with blocked polymerase activity at ... |
| <b>Form/Appearance</b> | liquid                                                                                                  |
| <b>Concentration</b>   | 2x conc                                                                                                 |
| <b>Storage</b>         | store at -20 °C                                                                                         |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                          |