


Product Datasheet

EvaGreen Fluorescent DNA Stain protein (orb533389)

Description	EvaGreen Fluorescent DNA Stain is a superior DNA intercalator dye specially developed for DNA analys...
Form/Appearance	liquid, supplied in 20 mM Tris-HCl pH 8.5, 0.1 mM EDTA and 0.01 % Tween-20; Color: orange
Concentration	100 µM
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{exc} 500 nm, λ_{em} 530 nm (bound to DNA).
Expiration Date	6 months from date of receipt.



Excitation
(black) and
emission
(green) ...