

## Product Datasheet

### Agarose Gel Extraction Kit - Column Kit (orb533437)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Agarose Gel Extraction Kit is designed for high-yield recovery of DNA from agarose gel with simultan... |
| <b>Storage</b>         | store at ambient temperature                                                                            |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                         |