

Product Datasheet

Atto425 NT Labeling Kit (orb533470)

Description	Atto425 Nick Translation Labeling Kit contains all reagents (except template and materials for purif...
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{exc} 436 nm, λ_{em} 484 nm, ϵ 45.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5).
SMILES	<chem>O=C(OCC)c1cc2c(oc1=O)cc1N(CCCC(=O)NC/C=C/c3cn(c(=O)[nH]c3=O)[C@H]3C[C@@H]([C@H](O3)COP(=O)(O)OP(=O)(O)OP(=O)(O)O)C(C)(C)CC(C)c1c2</chem>
Expiration Date	6 months from date of receipt.