

Product Datasheet

REP-1 His protein (orb533495)

Description	C-terminal His-tagged REP-1 (Rab escort protein-1) binds to a wide range of RabGGTase with nanomolar...
Form/Appearance	liquid (Supplied in 25 mM HEPES pH 7.2, 40 mM NaCl and 5 mM DTT)
Storage	store at -80 °C
Note	For research use only
Application notes	Activity: 1 pmol of REP-2 will mediate transfer of 2 pmol of geranylgeraniol by RabGGTase onto 1 pmol of Rab7 in 3 minutes.
Clonality	Recombinant
Purity	> 90% (SDS-PAGE)
Source	E. coli
MW	Theoretical MW: 75 kDa
Expiration Date	6 months from date of receipt.