

## Product Datasheet

### RabGDI protein (orb533500)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Rab GDP Dissociation Inhibitor (RabGDI) forms complexes with cytosolic prenylated GDP-bound Rab prot... |
| <b>Form/Appearance</b> | liquid (Supplied in 40 mM HEPES pH 7.2, 50 mM NaCl, 100 $\mu$ M GDP and 3 mM DTE)                       |
| <b>Storage</b>         | store at -80 °C                                                                                         |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Clonality</b>       | Recombinant                                                                                             |
| <b>Purity</b>          | > 90% (SDS-PAGE)                                                                                        |
| <b>Source</b>          | Sf9 insect cells                                                                                        |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                          |