

Product Datasheet

HAV-P2C (residues 1121-1234) protein (orb533503)

Description	The E. coli derived recombinant 38.5 kDa protein contains the P2C immunodominant regions, amino acids ...
Form/Appearance	liquid (Supplied in 10 mM CBB pH 9.6, 0.1% SDS and 50% glycerol)
Storage	store at -20 °C
Note	For research use only
Application notes	Recombinant HAV-P2C may be used in ELISA and Western blots, excellent for detection of HAV with minimal specificity problems.
Clonality	Recombinant
Purity	> 90% (SDS-PAGE)
MW	Theoretical MW: 38.5 kDa
Expiration Date	6 months from date of receipt.