

Product Datasheet

Rab2GST-His protein (orb533505)

Description	Rab proteins are small-molecular-weight guanosine triphosphatases (GTPases) that control vesicular t...
Form/Appearance	liquid (Supplied in 25 mM HEPES pH 7.2, 40 mM NaCl, 5 mM DTT, 2 mM MgCl 2and 10 mM GDP)
Storage	store at -80 °C
Note	For research use only
Clonality	Recombinant
Purity	> 90% (SDS-PAGE)
Source	E. coli
MW	Theoretical MW: 53.5 kDa
Expiration Date	6 months from date of receipt.