

Product Datasheet

HAV-VP4-VP2 (residues 55-164) protein (orb533506)

Description	The protein contains the HAV structural proteins VP2-VP4 immunodominant regions. Hepatitis A Virus ...
Form/Appearance	liquid (Supplied in 10 mM CBB pH 9.6, 0.1% SDS and 50% glycerol)
Storage	store at -20 °C
Note	For research use only
Application notes	Recombinant HAV-VP4-VP2 may be used in ELISA and Western blots, excellent for detection of HAV with minimal specificity problems.
Clonality	Recombinant
Purity	> 90% (SDS-PAGE)
MW	Theoretical MW: 44 kDa
Expiration Date	6 months from date of receipt.