

Product Datasheet

Rab17GST-His protein (orb533517)

Description	Rab17 is a small GTPase that belongs to the Ras superfamily. Rab proteins play an important role in ...
Form/Appearance	liquid (Supplied in 25 mM HEPES pH 7.2, 40 mM NaCl, 2 mM MgCl 2, 10 µM GDP and 3 mM DTE)
Storage	store at -80 °C
Note	For research use only
Clonality	Recombinant
Purity	> 90% (SDS-PAGE)
Source	E. coli
MW	Theoretical MW: 52 kDa
Expiration Date	6 months from date of receipt.