

Product Datasheet

Hepatitis B Virus Surface Antigen - adr subtype protein (orb533519)

Description	Recombinant HBsAg-adr full length monomeric protein contains 227 amino acids of the S-gene and has a...
Form/Appearance	liquid (Supplied in 20 mM sodium phosphate buffer (PB), 154 mM NaCl, pH 7.1)
Storage	store at 4 °C
Note	For research use only
Application notes	Antigen in ELISA, Immunochromatography and monoclonal antibody production.
Clonality	Recombinant
Purity	> 95% (RP-HPLC)
MW	Theoretical MW: 23-27 kDa
Expiration Date	6 months from date of receipt.