

Product Datasheet

West Nile Virus EnvelopeHis protein (orb533527)

Description	The E.coli derived recombinant protein contains the West-Nile N-terminus Envelope Virus immunodomain...
Form/Appearance	liquid (Supplied in 20 mM Phosphate buffer pH 7.5)
Storage	store at -20 °C
Note	For research use only
Application notes	Antigen in ELISA and Western Blots, excellent antigen for detection of West-Nile virus with minimal specificity problems.
Clonality	Recombinant
Purity	> 95% (SDS-PAGE)
Expiration Date	6 months from date of receipt.