

## Product Datasheet

### Parvovirus B19 VLP VP1/VP2 Co-Capsid protein (orb533528)

|                          |                                                                                                         |
|--------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Recombinant Parvovirus B19 VLP VP1/VP2 Co-Capsid is produced in SF9 containing B19 capsid proteins V... |
| <b>Form/Appearance</b>   | liquid (Supplied in 16 mM sodium phosphate pH 7.4, 120 mM NaCl and 20% glycerol)                        |
| <b>Storage</b>           | store at -20 °C                                                                                         |
| <b>Note</b>              | For research use only                                                                                   |
| <b>Application notes</b> | Western Blot with anti-Parvovirus B19 positive sera. Binds IgG-type human antibodies.                   |
| <b>Clonality</b>         | Recombinant                                                                                             |
| <b>Purity</b>            | > 80% (SDS-PAGE)                                                                                        |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                                          |