

Product Datasheet

Parvovirus B19 VLP VP2 protein (orb533529)

Description	Recombinant Parvovirus B19 VLP VP2 produced in SF9 is a glyco-sylated, polypeptide chain having a ca...
Form/Appearance	liquid (Supplied in 16 mM sodium phosphate pH 7.4, 120 mM NaCl and 20% glycerol)
Storage	store at -20 °C
Note	For research use only
Application notes	Western Blot with anti-Parvovirus B19 positive sera. Binds IgG-type human antibodies.
Clonality	Recombinant
Purity	> 80% (SDS-PAGE)
Expiration Date	6 months from date of receipt.