

Product Datasheet

Hantavirus GST protein (orb533530)

Description	The Hantavirus nucleocapsid (N) fusion protein is expressed in E.coli and fused to a GST-tag at the ...
Form/Appearance	liquid (supplied in 1 x PBS pH 7.4 and 0.05% sodium azide)
Storage	store at -20 °C
Note	For research use only
Application notes	Immunoassays to test the specific IgM and IgG to Hantavirus.
Clonality	Recombinant
Purity	> 95% (SDS-PAGE, 12%, coomassie staining)
Expiration Date	6 months from date of receipt.