

Product Datasheet

Chikungunya Wild Type E1 protein (orb533535)

Description	Recombinant Chikungunya Wild Type E1 produced in Insect Cells is a polypeptide chain containing amin...
Form/Appearance	liquid (Supplied in 1x PBS, pH 7.4, 0.1% Thimerosal, 5 mM EDTA, 1 µg/ml of Leupeptin, Aprotinin and Pepstatin A)
Storage	store at -20 °C
Note	For research use only
Clonality	Recombinant
Purity	> 95% (SDS-PAGE)
Expiration Date	6 months from date of receipt.