

## Product Datasheet

### Chikungunya Mutant (A226V) E1 protein (orb533536)

|                        |                                                                                                                  |
|------------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Recombinant Chikungunya Mutant (A226V) E1 produced in Insect Cells is a polypeptide chain containing...          |
| <b>Form/Appearance</b> | liquid (Supplied in 1x PBS, pH 7.4, 0.1% Thimerosal, 5 mM EDTA, 1 µg/ml of Leupeptin, Aprotinin and Pepstatin A) |
| <b>Storage</b>         | store at -20 °C                                                                                                  |
| <b>Note</b>            | For research use only                                                                                            |
| <b>Clonality</b>       | Recombinant                                                                                                      |
| <b>Purity</b>          | > 95% (SDS-PAGE)                                                                                                 |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                                   |