



Product Datasheet

TAF1A antibody (orb53355)

Description	TAF1A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human TATA box-binding protein-associated factor RNA polymerase I subunit A protein (251-325AA)
Target	TAF1A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q15573
Dilution Range	WB:1:200-3000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain tissue.



Western blot analysis of A549 whole cell lysate.