

Product Datasheet

Ran - T24N protein (orb533566)

Description	Ran (Ras-related nuclear protein, TC4) belongs to the superfamily of GTP binding proteins. It acts i...
Form/Appearance	liquid (Supplied in 20 mM KH ₂ PO ₄ /K ₂ HPO ₄ 4pH 7.0, 6.5 mM DTE, 150 mM KCl, 5 mM MgCl ₂ and 10% glycerol)
Storage	store at -80 °C
Note	For research use only
Clonality	Recombinant
Purity	> 90% (SDS-PAGE)
Source	E. coli
MW	Theoretical MW: 24.0 kDa (216 amino acids)
Expiration Date	6 months from date of receipt.