

Product Datasheet

Rac1 C protein (orb533576)

Description	Rho family GTPases Rac1 and Cdc42 belong to the Ras superfamily of small GTP-binding proteins. Rac1 ...
Form/Appearance	liquid (Supplied in 30 mM Tris-HCl pH 7.0 and 5 mM DTE)
Storage	store at -80 °C
Note	For research use only
Clonality	Recombinant
Purity	> 90% (SDS-PAGE)
Source	E. coli
MW	Theoretical MW: 20.9 kDa (187 amino acids)
Expiration Date	6 months from date of receipt.