

Product Datasheet

RafRBD (c-Raf) protein (orb533577)

Description	Raf is a member of the serine/threonine proteine kinase family involved in regulation of cell growth...
Form/Appearance	liquid (Supplied in 32 mM Tris-HCl pH 7.6, 200 mM NaCl, 5 mM MgCl 2and 5 mM DTT)
Storage	store at -80 °C
Note	For research use only
Clonality	Recombinant
Purity	> 90% (SDS-PAGE)
Source	E. coli
MW	Theoretical MW: 9.4 kDa (50-132 amino acids)
Expiration Date	6 months from date of receipt.