

## Product Datasheet

### RCC1 (RanGEF) protein (orb533578)

|                        |                                                                                                              |
|------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | RCC1 (Regulator of Chromosome Condensation, a $\beta$ propeller chromatin-bound protein) catalyzes guanin... |
| <b>Form/Appearance</b> | liquid (Supplied in 20 mM HEPES pH 7.6, 50 mM NaCl and 2 mM DTE)                                             |
| <b>Storage</b>         | store at -80 °C                                                                                              |
| <b>Note</b>            | For research use only                                                                                        |
| <b>Clonality</b>       | Recombinant                                                                                                  |
| <b>Purity</b>          | > 90% (SDS-PAGE)                                                                                             |
| <b>Source</b>          | E. coli                                                                                                      |
| <b>MW</b>              | Theoretical MW: 44.97 kDa (421 amino acids)                                                                  |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                               |