

Product Datasheet

FSME (TBE) Virus Antigen (orb533701)

Description	Strain Moscow B-4
Preservatives	Preservation: None
Form/Appearance	liquid
Storage	store at -80 °C
Note	For research use only
Application notes	<p>Culture system: Vero cell culture.</p> <p>Preparation: Supernatant from infected culture is concentrated by centrifugation. The resulting virus preparation is presented in glycine buffer at pH 9.5.</p> <p>Inactivation: This product has been treated with Beta-propiolactone. Confirmed inactive by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.</p>
Expiration Date	6 months from date of receipt.