

Product Datasheet

Brucella abortus Antigen (orb533704)

Description	S99 - Weybridge
Preservatives	Preservation: None
Form/Appearance	liquid (Supplied in glycine buffer pH 9.5 containing n-octyl- β -D-glucopyranoside)
Storage	store at -80 °C
Note	For research use only
Application notes	Culture system: Fermentation. Preparation: Bacterial lysate in glycine buffer (pH 9.5) contain-ing n-octyl-β-D-glucopyranoside. Inactivation: This product has been treated with heat (1 hour, 65°C), phenol (0.5% (w/w)) and detergent. Confirmed inactive by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.
Expiration Date	6 months from date of receipt.