

Product Datasheet

Helicobacter pylori Antigen (orb533706)

Description	Helicobacter pylori Antigen
Preservatives	Preservation: None
Form/Appearance	liquid
Storage	store at -80 °C
Note	For research use only
Application notes	Culture system: Solid medium. Preparation: Washed harvested organisms are solubilised and partially purified. Antigen is presented in a 1% non-ionic detergent solution. Inactivation: This product has been treated with detergent. Confirmed inactive by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.
Expiration Date	6 months from date of receipt.