

## Product Datasheet

### Parainfluenza Virus Type 1 Antigen (orb533715)

<b>Description</b>	Strain Sendai
<b>Preservatives</b>	Preservation: None
<b>Form/Appearance</b>	liquid
<b>Storage</b>	store at -80 °C
<b>Note</b>	For research use only
<b>Application notes</b>	<p><b>Culture system:</b> MDCK cell culture.</p> <p><b>Preparation:</b> Infected cell culture supernatant is treated to remove cellular material, further purified by size exclusion chro-matography and concentrated. The antigen is presented in Tris/NaCl buffer at pH 7.4.</p> <p><b>Inactivation:</b> This product has been treated with Beta-propiolactone. Confirmed inactive by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.</p>
<b>Expiration Date</b>	6 months from date of receipt.