

Product Datasheet

Parainfluenza Virus Type 3 Antigen (orb533717)

Description	Strain C243
Preservatives	Preservation: None
Form/Appearance	liquid
Storage	store at -80 °C
Note	For research use only
Application notes	<p>Culture system: Vero cell culture.</p> <p>Preparation: Infected cell culture supernatant is treated to remove cellular material, further purified by size exclusion chro-matography and concentrated. The antigen is presented in Tris/NaCl buffer at pH 7.4.</p> <p>Inactivation: This product has been treated with Beta-propiolactone. Confirmed inactive by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.</p>
Expiration Date	6 months from date of receipt.