

Product Datasheet

Yersinia enterocolitica O:9 Antigen (orb533727)

Description	Yersinia enterocolitica O:9 Antigen
Preservatives	Preservation: None
Form/Appearance	liquid (Supplied in PBS)
Storage	store at -80 °C
Note	For research use only
Application notes	<p>Culture system: Liquid culture.</p> <p>Preparation: Yersinia is grown in liquid culture and stimulated to secrete YOPs (Yersinia outer proteins). YOP proteins concentrated from the culture medium are presented in phosphate buffered saline.</p> <p>Inactivation: This product has been inactivated by sterile filtration. Confirmed inactive by attempted growth under original culture conditions. Since no procedure can guarantee absolute sterility, the reagent should be handled with appropriate safety precautions.</p>
Expiration Date	6 months from date of receipt.