

Product Datasheet

Amphipol A8-35 (orb533771)

Description	Amphipol A8-35 is a short amphipatic polymer that is specifically designed for membrane protein stab...
Storage	store at -20 °C
Note	For research use only
Application notes	Stabilizing agent in Cryo-EM and X-ray crystallography.
Formula	$(C_6 \cdot 2H_{10} \cdot 3O_1 \cdot 35N_0 \cdot 65Na_0 \cdot 35)_{72}$
MW	Theoretical MW: approx. 9 kDa
CAS Number	326856-53-5
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.