

Product Datasheet

Biorbyt Thermofluor Dye (orb533793)

Description	Biorbyt Thermofluor Dye
Form/Appearance	liquid
Concentration	5000x in DMSO
Storage	store at -20 °C
Note	For research use only
Application notes	Biorbyt Thermofluor Dye complements BiorbytScreen Thermofluor SPECIFIC HTS (CSS-333). Fluorescence readout is recommended in the range of 570-600 nm.
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.