

Product Datasheet

CD227 antibody (FITC) (orb533825)

Description	Mouse monoclonal antibody against to CD227 conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.09% sodium azide and 1% rAlbumin
Concentration	0.1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG3
Clonality	Monoclonal
Dilution Range	,
Expiration Date	12 months from date of receipt.