

## Product Datasheet

### Cadherin 16 Antibody (orb534027)

## Description

The CDH16 gene is a member of the cadherin

dependent, membrane-a...

## Species/Host

Rabbit

## Reactivity

Human, Mouse

## Conjugation

Unconjugated

## Tested

IHC-P, WB

## Applications

## Immunogen

Human protein was used as the immunogen for the recombinant Cadherin 16 antibody.

## Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

## Storage

Store the recombinant Cadherin 16 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

## Note

For research use only

## Application notes

The stated application concentrations are suggested starting points. Titration of the recombinant Cadherin 16 antibody may be required due to differences in protocols and secondary/substrate sensitivity.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

## Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

## Isotype

Rabbit IgG, kappa

## Clonality

Recombinant

## Purity

Protein A affinity chromatography

## Uniprot ID

**075309**

## Hazard

## Information

This recombinant Cadherin 16 antibody is available for research use only.

## Dilution Range

Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT, Prediluted IHC only format: incubate for 30 min at RT (1)

## Expiration Date

12 months from date of receipt.