

Product Datasheet

CD9 Antibody (orb534092)

Description

Involved in platelet activation and aggregation.

in...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

Recombinant full-length human protein was used as the immunogen for the CD9 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the CD9 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

The stated application concentrations are suggested starting points. Titration of the CD9 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P21926

Hazard
Information

This CD9 antibody is available for research use only.

Dilution Range

IHC (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.