

Product Datasheet

CD79b Antibody (orb534108)

Description

Required in cooperation with CD79A for

	activated by th...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A portion of amino acids 29-159 was used as the immunogen for the CD79b antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide
Storage	Store the CD79b antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	The stated application concentrations are suggested starting points. Titration of the CD79b antibody may be required due to differences in protocols and secondary/substrate sensitivity.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P40259
Hazard Information	This CD79b antibody is available for research use only.
Dilution Range	IHC (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.