

Product Datasheet

OLIG2 Antibody (orb534120)

Description

Olig2, a basic helix-loop-helix transcription

specification....

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

A portion of amino acids 1-141 from the human protein was used as the immunogen for this OLIG2 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

After reconstitution, the OLIG2 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Optimal dilution of the OLIG2 antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity

Uniprot ID

Q13516

Hazard

Information

This OLIG2 antibody is available for research use only.

Dilution Range

IHC (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.